



Release Notes

ARUM 3D Dental Solutions library for 3Shape

Version 1.2.0 [21/03/2025]

[New Compatible System Library Update]

Pre-mill & Disk

MEDENTiKA® Multi-Unit

[Issue Fixes and Updates]

Pre-mill & Disk/Ti-base

ASTRA TECH™ EV™: Corrected 60-degree Atlantis scan body rotation alignment

OSSTEM® SS Regular 4.8: Adjusted Arum scan body height position

Pre-mill & Disk

EUROTEKNIKA™ Multi-Unit: Updated EURO-ETK-UN-SBP scan body file

Version 1.1.2 [03/03/2025]

[New Compatible System Library Update]

Pre-mill & Disk

Global D® EVL® Multi-Unit ST

Global D® EVL® Multi-Unit WD

[Ti-Base Library Update]

BTI® INTERNA®



KENTEC®

NEOSS Pro Active®

Version 1.0.2 [07/02/2025]

[Issue Fixes and Updates]

We have identified an issue in the recently distributed 3Shape CAD library where the **D5.0 and D6.0 analog files** have size discrepancies, causing improper fitting during printed model design.